

PANDA: Portable and Adaptable Neutron Diagnostics for ARPCA-E

**FUSION Program Review (Virtual)
March 5, 2021**

Drew P. Higginson, LLNL

LLNL: Chris Cooper, Clement Goyon, Matt McMahon, James Mitrani,
Amanda Youmans



UCB/LBNL: Bethany Goldblum (PI), Josh Brown, Thibault Laplace



U.S. DEPARTMENT OF
ENERGY

Team members and roles

LLNL



UCB/LBNL


LLNL Team

- **Drew Higginson** – *PI*
- **James Mitrani** – *Detector design, assembly and calibration*
- **Amanda Youmans** – *Detector assembly, calibration and fielding, data analysis*
- **Clement Goyon** – *LaBr analysis*
- **Chris Cooper** – *LaBr design and assembly*
- **Matt McMahon** – *MCNP analysis of neutron scattering*

UCB/LBNL Team

- **Bethany Goldblum** – *PI*
- **Josh Brown** – *Geant4 simulations and data analysis*
- **Thibault Laplace** – *Pulse metrology and data analysis*

Our diagnostics measure neutron yield and infer neutron energy for neutron emission as small as a few ns to multi-millisecond.





Total neutron yield via activation:


- Viable at yields >5e6 neutrons.
- Yield calculated to better than 20%.

Thermonuclear & beam-target fusion discrimination via Scintillator/PMTs:

- Viable at yields >1e5.
- Few ns emission: time-of-flight method.
- Multi- μ s emission: pulse counting method
- Ion beams > 100 keV can be detected.




Diagnostics ship in single EMI racks that are easy to transport.




- Simple, compact, easy to transport & non-invasive.
- Detectors: **3x** LaBr activation & **24x** S/PMTs.
- Easy to use scripts; provide raw and analyzed data for the user. Data delivered <3 minutes from a shot.
- We can field and analyze data. We will also teach the teams to run themselves.

LaBr Analysis Script Outputs




Activation & S/PMTs calibrated & tested on the LLNL Mjolnir DPF

Activation Detector Calibration Setup




Activation Calibration Results




Distance	Min Yield
20 cm	5.0e+06
50 cm	3.1e+07
100 cm	1.2e+08

S/PMT Test Setup



- Calibration and test shots taken on the Mjolnir DPF at LLNL
- Diagnostics performed very well in this harsh environment.
- All systems ready for deployment.




Ready to deploy! MIFTI (UC San Diego) & ZEI (Zap lab)


MIFTI at UC San Diego

- MCNP sims of device/room used to determine the best diagnostic placement .
- Shipping in March; deuterium ops coming soon.

H-factor = expected neutrons / neutrons in vacuum




Close-up on device




Zap Energy Inc. at Zap Lab

- Monte-Carlo: ions \rightarrow neutrons \rightarrow protons \rightarrow pulses
- Used data pipeline w/ previous FuZE data from ALPHA.
- Plan to take data at first plasma, \sim May 2021.

Geant4 Model



Data vs Model



Future Plans

- We aim to work closely with teams at MIFTI and ZEI to fully realize the data analysis methods and uncertainties. This will result in publications to help make the method clear to others.
- We will work with the teams if they want to bring up their own versions of these diagnostics (e.g., procurement, calibration, analysis scripts).
- This work has already been interesting to LLNL's DPF where anisotropy is expected and understanding ion beam dynamics is a core question for device optimization.